**DESCRIPTION:**

Nanohybrid resin for making dental crowns in resin.

**PRESENTATIONS:**

- Bottle containing 300 g of PriZma 3D Bio Crown Diamond resin.
  - Bottle containing 500 g of PriZma 3D Bio Crown Diamond resin.
- Available colors: BL, B1, B2, B3, A1, A2, A3, C1, C2, C3, PE, DBL, DA1, DA2, DA3, DB1, DB2, DB3, PEBL.

**INDICATION OF USE:**

PriZma 3D Bio Crown Diamond is an acrylic resin suitable for 3D printing to make resin dental crowns.

**COMPOSITION:**

Oligomers, monomers, photoinitiators, UV absorbers, fillers, and pigments.

**MODE OF USE:**

1. Wear gloves and goggles to handle the product.
2. Put the amount of resin indicated by the printer manufacturer in the printer tank and start printing.
3. After printing, put the objects in isopropyl alcohol for 3 to 5 minutes in motion to remove the resin from the surface and dry completely.
4. Put in UV light chamber for post-curing for 5 to 25 minutes. Polish/glaze and wash before taking to mouth.
5. Do not leave the resin in the tank after printing.
6. Optimal working temperature on the printer from 25 to 38 °C. Perform resin calibration at each batch change.
7. Check compatible printers (3D printers with wavelengths between 295 to 410 nm).

**WARNINGS:**

- The product must be handled by a specialized professional.
- Do not use if the expiration date has passed.
- The product must be handled with gloves, goggles and, preferably, in a ventilated environment.
- The product can be irritating to the eyes, skin, and respiratory system.
- Keep the material out of the reach of children.
- This material is manufactured for dental use only and must be handled in accordance with the instructions for use. The manufacturer is not responsible for damage caused by other uses or by incorrect handling not described in this manual. In addition, the user is obliged to check, before use and under their responsibility, whether this material is compatible with the intended use, especially when this use is not indicated in this indication of use.
- Strictly follow the procedures provided in the instructions for use.
- Do not use contaminated spatulas to remove material from the inside of the 3D printer.
- Contact of the non-light-cured PriZma 3D Bio Crown Diamond with the skin, mucous membranes and eyes should be avoided. When not yet polymerized, the product may cause a slightly irritating effect.
- In case of allergic reactions, discontinue use immediately. If necessary, seek medical advice.
- Single-use product. After photoactivation, the product does not return to its initial consistency.
- Reprocessing is prohibited.
- Do not expose the material directly to sunlight.
- After use, the perfect sealing of the bottle must be ensured to prevent the material from being exposed to light.
- Avoid contact with the eyes so as not to damage the cornea. In case of reaction with this fabric, you should wash with plenty of water and seek medical assistance.
- Avoid contact with the oral mucosa in order to avoid the formation of inflammatory processes. In accidental contacts, washing with plenty of water and medical assistance is recommended.
- Avoid ingesting the product. In accidental cases, the patient should drink plenty of water and, in case of unwanted reactions, medical assistance should be

sought.

- Do not use the material if the packaging is damaged.
- Do not perform post-curing for longer than recommended by the manufacturer.

**PRECAUTIONS:**

- Avoid contamination of land, plumbing and surface water.
- Toxic to aquatic organisms, it can cause long-term adverse side effects on the aquatic environment.
- Keep out of reach of children and animals.
- This product should only be used in accordance with the instructions for use.
- Wear PPE when handling the product, including gloves and goggles. It is recommended that patients also wear protective goggles.
- In case of storage in refrigerated systems, before application, the material should be left at room temperature so that it reaches adequate fluidity.

**FIRST AID**

**Contact with skin:** Wipe with absorbent paper or towel. Wash immediately with plenty of soap and water. Do not use organic solvents. In the case of dermatitis, seek medical help.

**Contact with eyes:** Rinse with plenty of running water for 20 to 30 minutes. Refer to doctor if necessary.

**Inhalation:** Go to a cool and ventilated place if there is strong inhalation of the product, call/refer to the doctor if necessary.

**Ingestion:** Do not induce vomiting. The affected person should drink 0.5 – 0.8 L of water and seek medical help.

For more information, see the SDS available on the [www.makertechlabs.com.br](http://www.makertechlabs.com.br).

**ADVERSE EFFECTS:**

So far, no adverse effects of the product have been identified. No side effects are expected if the product is used as directed.

**STORAGE CONDITIONS:**

Store PriZma 3D Bio Crown Diamond resin in its original packaging and in a clean, dry place. Keep the material properly sealed and without exposure to light or high humidity. The product should be kept at a controlled temperature between 5 °C and 30 °C.

**TRANSPORT CONDITIONS:**

PriZma 3D Bio Crown Diamond resin can be transported for 15 days between extreme temperatures of -28 °C and 55 °C without compromising its integrity. Do not transport at temperatures outside the recommended range. Do not place heavy objects on the product packaging.

**VALIDITY:**

If stored as instructed by the manufacturer, the product has a shelf life of 24 months from the date of manufacture.

Once opened, if stored according to the manufacturer's recommendations, the product remains stable for 3 months.

**DISPOSAL:**

Dispose of the product in accordance with local legislation.

It is recommended that the professional apply good clinical practices for the disposal of medical products according to current health standards.

Do not use after the expiration date.

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