DESCRIPTION:

A UVA/UVB-curable acrylic resin designed exclusively for use in 3D printers, suitable for producing myorelaxant splints with enhanced toughness to ensure greater comfort and durability.

PRESENTATIONS:

- Bottle with 250 g of PriZma 3D Bio Splint Soft resin;
- Bottle with 500 g of PriZma 3D Bio Splint Soft resin;
- Bottle with 1000 g of PriZma 3D Bio Splint Soft resin.

INDICATION FOR USE:

PriZma 3D Bio Splint Soft resin is indicated exclusively for use in 3D printers for manufacturing myorelaxant plaques. Myorelaxant splints are intended to protect teeth from wearing out in patients with bruxism (habit of grinding teeth), relieve temporomandibular joints and induce muscle relaxation.

COMPOSITION:

Oligomers, Monomers and Pigments.

MODE OF USE:**1. Printing**

- Add the amount of resin indicated by the manufacturer to the printer tank.
- Wear gloves and goggles throughout handling.
- Start printing according to the printer parameters.

2. Washing

- After printing, wash the parts for 3 to 5 minutes on the go, to remove surface resin.
- Use water with dental detergent/descaler and/or 70% alcohol.
- Make sure the piece is completely dry before the next step.

3. Finishing and Pre-Curing

- After drying, remove the supports and finish the surface with finishing cutters until it is smooth.
- Clean the piece and apply a thin layer of the resin itself to the entire surface with a brush.
- Perform a 1-minute pre-cure in a UV light chamber.

Attention: avoid excess resin on the inside so as not to compromise the adaptation to the teeth.

4. Post-Curing & Shine

- Submerge the piece in glycerin and perform post-curing for 8 to 15 minutes in UV light. This process maximizes the natural shine of the piece, eliminating mechanical polishing. If you opt for mechanical polishing, the result may be limited due to the soft characteristic of the resin.

Important Notes:

- Never leave resin stored in the tank after printing.
- Optimal printer operating conditions: 25 °C to 38 °C.

WARNINGS:

- The product must be handled by a specialized professional.
- Do not use if the expiration date has expired.
- The product must be handled with gloves, goggles and, preferably, in a ventilated environment.
- The product can be irritating to the eyes, skin, and respiratory system.
- Keep the material out of the reach of children.
- This material has been manufactured for dental use and must be handled in accordance with the instructions for use. The manufacturer is not responsible for damage caused by other uses or by incorrect handling not described in this manual. In addition, the user is obliged to check, before use and under his responsibility, whether this material is compatible with the intended use, especially when this use is not indicated in this indication of use.
- Strictly follow the procedures provided in the instructions for use.
- Do not use contaminated spatulas to remove material from the inside of the 3D printer.
- Contact of the non-light-cured PriZma 3D Bio Splint Soft with the skin, mucous membranes, and eyes should be avoided. When not yet polymerized, the product may cause a slightly irritating effect.
- In case of allergic reactions, discontinue use immediately. If necessary, seek medical advice.
- Single-use product. After photoactivation, the product does not return to its initial consistency.
- Reprocessing is prohibited.
- Do not expose the material directly to sunlight.

- After use, the perfect sealing of the bottle must be observed to prevent the material from being exposed to light.
- Avoid contact with the eyes so as not to damage the cornea. In case of reaction with this tissue, you should wash with plenty of water and seek medical assistance.
- Avoid contact with the oral mucosa in order to avoid the formation of inflammatory processes.
- In accidental contacts, washing with plenty of water and medical assistance is recommended.
- Avoid ingestion of the product. In accidental cases, the patient should drink plenty of water and, in case of unwanted reactions, medical assistance should be sought.
- Do not use the material if the packaging is damaged.
- Do not perform post-curing for longer than recommended by the manufacturer.
- Do not leave submerged for long periods in the liquids indicated for washing.
- It is not recommended to use any other product for cleaning or submerging the piece in any liquid not indicated in this manual.

PRECAUTIONS:

- This product should only be used in accordance with the instructions for use.
- Wear PPE when handling the product, including gloves and goggles. It is recommended that patients also wear protective goggles.
- In case of storage in refrigerated systems, before application, the material should be left at room temperature so that it reaches adequate fluidity.

FIRST AID:

CONTACT WITH SKIN: Wipe with absorbent paper or towel. Wash immediately with plenty of soap and water. Do not use organic solvents. In the case of dermatitis, seek medical help.

CONTACT WITH EYES: Rinse with plenty of running water for 20 to 30 minutes. Refer to doctor if necessary.

INHALATION: Go to a cool and ventilated place if there is strong inhalation of the product, call/refer to the doctor if necessary.

INGESTION: Do not induce vomiting. The affected person should drink 0.5 – 0.8 L of water and seek medical help.

For more information, see the SDS available on www.makertechlabs.com.br.

CONTRAINDICATION:

PriZma 3D Bio Splint Soft should not be used in case of known hypersensitivity (allergy) to any of the components of the formula.

ADVERSE EFFECT:

So far, no adverse effects of the product have been identified. However, there are reports in the literature of adverse effects with monomers and oligomers. It is the dentist's responsibility to check the allergic condition of each patient.

STORAGE CONDITIONS:

Store PriZma 3D Bio Splint Soft resin in its original packaging and in a clean, dry place. Keep the material properly sealed and without exposure to light or high humidity. The product should be kept at a controlled temperature between 5 °C and 30 °C.

TRANSPORT CONDITIONS:

PriZma 3D Bio Splint Soft resin can be transported for 15 days in extreme temperatures of -28 °C and 55 °C without compromising its integrity. Do not transport at temperatures outside the recommended range. Do not place heavy objects on the product packaging.

VALIDITY:

If stored as instructed by the manufacturer, the product has a shelf life of 24 months from the date of manufacture. Once opened, if stored according to the manufacturer's recommendations, the product remains stable for 3 months. Do not use after the expiration date.

MODE OF DISPOSAL:

Dispose of the product in accordance with local legislation. It is recommended that the professional apply good clinical practices for the disposal of medical products according to current health standards.

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