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DESCRIPTION:
PriZma 3D Makeup resin is a light-curing pigment available in different colors to allow makeup on composite resins and printed pieces.

PRESENTATIONS:
PriZma 3D Makeup available in colors: Intense A, Intense B, Terracotta, Honey, Vanilla, Purple, Acqua, Lunar, Red, Pink, Yellow and Black.

- 1 Syringe of PriZma 3D Makeup with 2.5 g.
- Kit with 8 PriZma 3D Makeup resin pigments with 2.5 g each and 1 bottle of PriZma 3D Glaze with 4 mL.

INDICATION:
PriZma 3D Makeup is indicated for:

- Makeup of provisionals printed in 3D resin;
- Color adaptation of teeth in full prosthetics;
- Color correction of artificial teeth;
- Characterization of partial dentures and artificial gums;
- Characterization of ceramic-zirconia reinforced composites, conventional composites and hybrid ceramics.

COMPOSITION:
Methacrylic monomers, initiators, co-initiators, preservatives, fillers and pigments.

MODE OF USE:
General Care
Before applying PriZma 3D Makeup light-curing pigments, the surface must be properly prepared and the occlusion and contact points must be adjusted.

- Application**
1. Aluminum oxide should be blasted onto the part to be characterized to increase the mechanical adhesion between the part and the pigment. In addition to blasting, proper cleaning should be carried out with isopropyl alcohol.
 2. Apply PriZma 3D Makeup resin pigment using a soft bristle brush until the desired effects and colors are achieved. PriZma 3D Glaze can be used as a thinning agent.
 3. Once the desired colors have been obtained, carry out the polymerization properly (in a UV chamber for 10 to 15 minutes or in a light-curing device for 8 minutes).
 4. After the pigment has polymerized, PriZma 3D Glaze should be applied. To do this, apply an even layer of PriZma 3D Glaze with a brush, in a single direction to avoid the incorporation of bubbles.
 5. Light-cure in UV light chamber for 10 to 15 minutes or light-curing appliance for 8 minutes (scanning every 4 sets of teeth).
 6. Apply the second layer and light-cure again under the same conditions as described above. It is advisable to apply a layer of oxygen inhibitor to ensure complete polymerization. After application, the material should be photoactivated again to increase the polymerization of the printed piece.

- WARNINGS**
- For the exclusive use of dentists and prosthetic technicians.
 - Do not reuse empty packaging.
 - Keep the material out of the reach of children.
 - This material has been manufactured for dental use only and must be handled in accordance with the instructions for use. The manufacturer is not responsible for damage caused by other uses or by incorrect handling not described in this manual. In addition, the user is obliged to check, prior to use and under their own responsibility, that this material is compatible with the intended use, especially when this use is not indicated in the present instructions for use.
 - Strictly follow the procedures provided in the instructions for use.
 - Do not allow loose polymerized objects and parts to be swallowed or aspirated.
 - Contact of the unpolymerized product with the skin, mucous membranes and eyes should be avoided. When not yet polymerized, it can have a slightly irritating effect.
 - In the event of allergic reactions, discontinue use immediately. If necessary, seek medical advice.
 - Single use product. After photoactivation, the product does not return to its original consistency.
 - Reprocessing forbidden.
 - Do not expose the material to temperatures above 25 °C.
 - Do not expose the material directly to sunlight.
 - After use, make sure the packaging is perfectly sealed to prevent the material from being exposed to light.
 - Avoid contact with the eyes to not damage the cornea. In the event of a reaction with this tissue, rinse with plenty of water and seek medical attention.
 - Avoid contact with the oral mucosa in order to prevent the formation of inflammatory processes. In case of accidental contact, rinse with plenty of water and seek medical attention.
 - Avoid ingesting the product. In accidental cases, the patient should drink plenty of water and, in the event of unwanted reactions, medical assistance should be sought.
 - Do not use the material if the packaging is damaged.

- PRECAUTIONS:**
- This product should only be used in accordance with the instructions for use.
 - Wear PPE when handling the product, including gloves and goggles. It is recommended that patients also wear goggles.
 - In the case of storage in refrigerated systems, the material should be left at room temperature before application so that it reaches the appropriate fluidity.
 - Light-cure in a UV light chamber for 10 to 15 minutes or in a light-curing appliance for 8 minutes.
 - Each UV chamber has a different intensity and a time of 10 to 15

minutes should be tested before applying the product. If this time is not satisfactory for the clinical case, it should be left in the UV chamber for longer to ensure complete polymerization.

- Light curing must be carried out properly in order to guarantee the product's durability and color stability.
- The professional must have experience in characterizing prosthetic parts and dental medical devices in order to ensure that the product achieves its indication and performance.
- Do not allow PriZma 3D Makeup resin pigment to come into contact with eugenol-based products, as polymerization inhibition may occur.

CONTRAINDICATION:
PriZma 3D Makeup should not be used in case of known hypersensitivity (allergy) to any of the components of the formula.

ADVERSE EFFECT:
So far, no adverse effects of PriZma 3D Makeup have been identified. However, there are reports in the literature of allergic reactions involving dental products based on methacrylic monomers.

STORAGE CONDITIONS:
Store PriZma 3D Makeup in its original packaging in a clean, dry place. Keep the material properly sealed and away from light or high humidity. The product should be kept at a controlled temperature between 4 °C and 25 °C.

TRANSPORT CONDITIONS:
The product can be transported between -28 °C and 55 °C for up to 15 days without compromising its initial performance.

VALIDITY:
If stored according to the manufacturer's instructions, the product has a shelf life of 24 months from the date of manufacture. Do not use after expiration date.

DISPOSAL:
Do not reuse the product and dispose of in accordance with local legislation.

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makertechLabs®
3D Solutions

prizma®
MAKEUP

 **MADE IN BRAZIL**  **maquira®**
Dental Group

 **Fabricante | Manufacturer:**
Maquira Indústria de Produtos Odontológicos S.A.
CNPJ: 05.823.205/0001-90
Endereço | Dirección | Address:
Av. Melvin Jones, 773
Pq. Ind. Bandeirantes
Maringá - PR | 87070-030 | Brasil
Resp. Técnico | Technical Resp.:
Cristiane Furlan | CRF-PR: 14.038
Reg. ANVISA Nº: 80322409026

 **SAP:**
Atendimento ao cliente
Atención al cliente
Customer service
0800 726 5848
sap@maquira.com.br