**DESCRIPTION:**

Nanohybrid resin for preparation of dental crowns.

**PRESENTATION:**

- Bottle with 300 g of Prizma 3D Bio Crown Diamond resin.
  - Bottle with 500 g of Prizma 3D Bio Crown Diamond resin.
- Available colors: BL, B1, B2, B3, A1, A2, A3, C1, C2, C3, PE, DBL, DA1, DA2, DA3, DB1, DB2, DB3 and PEBL.

**INDICATION OF USE:**

Prizma 3D Bio Crown Diamond is an acrylic resin suitable for 3D printing to make resin dental crowns.

**COMPOSITION:**

Monomers, photoinitiators, fillers and pigment.

**MODE OF USE:**

- Wear gloves and goggles when handling the product.
- Pour the amount indicated by the printer manufacturer into the printer tank and start printing.
- After printing, place the objects in isopropyl alcohol for 3 to 5 minutes while moving to remove the resin from the surface and dry completely.
- Place in a UV light chamber to post-cure for 5 to 25 minutes. Polish/glaze and wash before putting in mouth.
- Do not leave the resin in the tank after printing.
- Ideal working temperature in the printer is 25 to 38 °C. Calibrate the resin with each batch change.
- Check compatible printers (3D printers with a wavelength between 295 and 410 nm).
- Do not reuse the product.

**WARNINGS:**

- The product must be handled by a specialized professional.
- The product must be used in accordance with the instructions for use.
- The product must not be used if the expiry date has passed.
- The product must be handled with gloves, goggles and preferably in a ventilated environment.
- This product can be irritating to the eyes, skin and respiratory system.

**Skin contact:** Wipe off with absorbent paper or towel. Wash immediately with plenty of soap and water. Do not use organic solvents. In the event of dermatitis, seek medical attention.

**Eye contact:** Rinse with plenty of running water for 20 to 30 minutes. Refer to a doctor if necessary.

**Inhalation:** Move to fresh air in case of strong

inhalation of the product, call/refer a doctor if necessary.

**Ingestion:** Do not induce vomiting. The affected person should drink 0.5 - 0.8 L of water and seek medical attention. For further information, consult the SDS/INSTRUCTIONS available at [www.makertechlabs.com.br](http://www.makertechlabs.com.br).

**PRECAUTIONS:**

- Store in a dry, ventilated place, away from sources of heat and light.
- Keep container closed when not in use.
- Avoid contact with skin, eyes and clothing. May cause sensitization through skin contact.
- Do not ingest or inhale.
- Wear protective gloves and goggles when handling the product.
- Avoid contamination of land, drains and surface water.
- Toxic to aquatic organisms, may cause long-term adverse side effects in the aquatic environment.
- Keep out of reach of children and animals.

**ADVERSE EFFECT:**

To date, no adverse effects of the product have been identified. No side effects are expected if the product is used according to the instructions.

**STORAGE CONDITIONS:**

Store in a dry, ventilated place, away from sources of heat and light. Store the product in its original packaging, always closed when not in use. Store at temperatures between 10 and 35 °C.

**TRANSPORT CONDITIONS:**

Transport in its original packaging and at a temperature between 10 and 45 °C. Do not place heavy objects on the product packaging.

**VALIDITY:**

The product has a shelf life of 2 years after manufacture. If stored according to the manufacturer's instructions, the product can be used for up to 3 months after opening. Do not use after expiration date.

**DISPOSAL:**

Dispose of in accordance with local legislation.

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Nanohybrid resin for  
preparation of dental  
crowns

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