DESCRIPTION:

PriZma 3D Bio Guide is an acrylic resin indicated for printing/making surgical guides, orthognathic surgery guides and implant positioning guides.

PRESENTATIONS:

Bottle with 250 g of PriZma 3D Bio Guide resin.
Bottle with 500 g of PriZma 3D Bio Guide resin.
Bottle with 1000 g of PriZma 3D Bio Guide resin.

INDICATION FOR USE:

The PriZma 3D Bio Guide resin is a UVA/UVB light curable resin exclusively for use in 3D printers indicated for surgical guides for implants, gingival plasty, gutters for orthognathic surgeries and others.

COMPOSITION:

Monomers and photoinitiators.

MODE OF USE:

- Put the amount of resin indicated by the printer manufacturer in the printer tank and start printing;
- Wear gloves to handle the product;
- After printing, wash the objects for 3 to 5 minutes in motion to remove resin from the surface, wash with water and dental detergent/descaler and dry completely;
- Put into UV light chamber for post-curing: 30-60 minutes;
- Do not leave the resin in the tank after printing;
- Optimum working temperature on the printer from 25 °C to 38 °C.

NOTE: After handling and preparation, the surgical guide must be sterilized before use. Before use, the surgical guide must be placed in appropriate packaging for sterilization (surgical grade) and sterilized in an autoclave at 121 °C, according to the autoclave manufacturer's recommendations.

WARNINGS:

- The product must be handled by a specialized professional.
- Do not use if the expiration date has passed.
- The product must be handled with gloves, goggles and, preferably, in a ventilated environment.
- The product can be irritating to the eyes, skin, and respiratory system.
- Keep the material out of the reach of children.
- This material is manufactured for dental use only and must be handled in accordance with the instructions for use. The manufacturer is not responsible for damage caused by other uses or by incorrect handling not described in this manual. In addition, the user is obliged to check, before use and under his responsibility, whether this material is compatible with the intended use, especially when this use is not indicated in this indication of use.
- Strictly follow the procedures provided in the instructions for use.
- Do not use contaminated spatulas to remove material from the inside of the 3D printer.
- Contact of the non-light-cured PriZma 3D Bio Guide with the skin, mucous membranes and eyes should be avoided. When not yet polymerized, the product may cause a slightly irritating effect.
- In case of allergic reactions, discontinue use immediately. If necessary, seek medical advice.
- After photocuring, the product does not return to the initial consistency.
- Reprocessing is prohibited.
- Single-use product.
- Do not expose the material directly to sunlight.
- After use, the perfect sealing of the bottle must be observed to prevent the material from being exposed to light.
- Avoid contact with the eyes so as not to damage the cornea. In case of reaction with this tissue, you should wash with plenty of water and seek medical assistance.
- Avoid contact with the oral mucosa in order to avoid the formation of inflammatory

processes. In accidental contacts, washing with plenty of water and medical assistance is recommended.

- Avoid ingesting the product. In accidental cases, the patient should drink plenty of water and, in case of unwanted reactions, medical assistance should be sought.
 - Do not use the material if the packaging is damaged.
 - Do not perform post-curing for longer than recommended by the manufacturer.
- NOTE:** After handling and preparation, the surgical guide must be sterilized before use. Before use, the surgical guide must be placed in appropriate packaging for sterilization (surgical grade) and sterilized in an autoclave at 121 °C, according to the autoclave manufacturer's recommendations.

PRECAUTIONS:

- This product should only be used in accordance with the instructions for use.
- Wear PPE when handling the product, including gloves and goggles. It is recommended that patients also wear protective goggles.
- In case of storage in refrigerated systems, before application, the material should be left at room temperature so that it reaches adequate fluidity.

FIRST AID:

CONTACT WITH SKIN: Wipe with absorbent paper or towel. Wash immediately with plenty of soap and water. Do not use organic solvents. In the case of dermatitis, seek medical help.

CONTACT WITH EYES: Rinse with plenty of running water for 20 to 30 minutes. Refer to doctor if necessary.

INHALATION: Go to a cool and ventilated place if there is strong inhalation of the product, call/refer to the doctor if necessary.

INGESTION: Do not induce vomiting. The affected person should drink 0.5 – 0.8 L of water and seek medical help.

For more information, see the SDS available on the www.makertechlabs.com.br website.

CONTRAINDICATION:

PriZma 3D Bio Guide should not be used in case of known hypersensitivity (allergy) to any of the components of the formula.

ADVERSE EFFECT:

Further, no adverse effects of the product have been identified. However, there are reports in the literature of adverse effects with monomers and oligomers. It is the dentist's responsibility to check the allergic condition of each patient.

STORAGE CONDITIONS:

Store PriZma 3D Bio Guide resin in its original packaging and in a clean, dry place. Keep the material properly sealed and without exposure to light or high humidity. The product should be kept at a controlled temperature between 5 °C and 30 °C.

TRANSPORT CONDITIONS:

The product can be transported for 15 days in extreme temperatures between -28 °C and 55 °C without prejudice to its integrity. Do not transport at temperatures outside the recommended range. Do not place heavy objects on the product packaging.

VALIDITY:

If stored as instructed by the manufacturer, the product has a shelf life of 24 months from the date of manufacture. Once opened, if stored according to the manufacturer's instructions, the product can be used for up to 3 months. Do not use after the expiration date.

MODE OF DISPOSAL:

Dispose of the product in accordance with local legislation. It is recommended that the professional apply good clinical practices for the disposal of medical products according to current health standards.

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Resin for the manufacture
of surgical guides for
implants and others.

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