DESCRIPTION:
PriZma 3D Glaze is a light-curing varnish used on resins printed on 3D printers, recommended for sealing any porosity and giving the material the desired shine. This seal stabilises the surface, is resistant to abrasion and provides high shine.

PRESENTATION:
1 bottle of PriZma 3D Glaze with 4 or 6 mL.

INDICATION FOR USE:
Monomer-based varnish indicated for covering porosity in resins for printing total and partial prostheses on a 3D printer, sealing the surface and providing greater abrasion resistance and the desired shine to the material.

COMPOSITION:
Methacrylic monomers, initiators, co-initiators, stabilisers, blocker, fluorescent agent and filler.

MODE OF USE:
In the dental office
1. Prepare the material properly.
2. Apply a thin layer of the product using a soft brush. It is essential that the application be carried out in a single direction in order to avoid the formation of bubbles.
3. For each layer applied, light cure for 40 seconds using a scanning motion.

Resin thinner for characterisation
1. After curing the printed piece, sandblast it with aluminium oxide to increase adhesion.
2. Clean the part with isopropyl alcohol.
3. Dry the part thoroughly.
4. Apply thin layers of resin pigments to the surface to be characterised, moistening the brush with as little thinner as possible, just to facilitate application and allow the resin to be modelled.
5. Polymerise each layer and finish the characterisation according to the material specification.
6. After achieving the desired colour, apply a uniform layer of PriZma 3D Glaze with a brush, in a single direction to avoid incorporating bubbles.
7. Light-cure in a UV light chamber for 10 to 15 minutes or in a light-curing unit for 4 minutes (scanning every 4 sets of teeth).
8. Apply the second layer and light cure again under the same conditions as described above. It is advisable to apply a layer of oxygen inhibitor to ensure complete polymerization. After application, the material should be light cured again to increase the polymerisation of the printed piece.

WARNINGS:
- For exclusive use by dental surgeons and dental technicians.
- Do not reuse empty packaging.
- Keep the material out of the reach of children.
- This material has been manufactured for dental use only and must be handled in accordance with the instructions for use. The manufacturer is not liable for damage caused by other uses or by incorrect handling not described in this manual. In addition, the user is obliged to verify, prior to use and under their own responsibility, whether this material is compatible with the desired use, especially when this use is not indicated in the present instructions for use.
- Strictly follow the procedures provided in the instructions for use.
- Do not allow loose polymerized objects and parts to be swallowed or inhaled.
- Avoid contact of the unpolymerized product with the skin, mucous membranes and eyes. When not yet polymerized, it may cause a slightly irritating effect.
- In case of allergic reactions, discontinue use immediately. If necessary, seek medical advice.
- Single-use product. After photoactivation, the product does not return to its initial consistency.
- Reprocessing is prohibited.
- Do not expose the material to direct sunlight.
- After use, ensure that the packaging is completely sealed to prevent the material from being exposed to light.
- Avoid contact with the eyes to prevent damage to the cornea. In case of reaction with this tissue, wash with plenty of water and seek medical assistance.
- Avoid contact with the oral mucosa to prevent the formation of inflammatory processes. In case of accidental contact, rinse with plenty of water and seek medical assistance.
- Avoid ingesting the product. In case of accidental ingestion, the patient should drink plenty of water and, in case of adverse reactions, seek medical assistance.
- Do not use the material if the packaging is damaged.

PRECAUTIONS:
- This product should only be used in accordance with the instructions for use.
- Wear PPE when handling the product, including gloves and safety glasses. It is recommended that patients also wear safety glasses.
- If stored in refrigerated systems, the material should be left at room temperature before application so that it reaches the appropriate fluidity.

CONTRAINDICATIONS:
PriZma 3D Glaze should not be used in cases of known hypersensitivity (allergy) to any of the components of the formula.

ADVERSE EFFECTS:
To date, no adverse effects of PriZma 3D Glaze have been identified. However, there are reports in the literature of allergic reactions involving dental products based on methacrylic monomers.

STORAGE CONDITIONS:
Store PriZma 3D Glaze in its original packaging in a clean, dry place. Keep the material properly sealed and away from light or high humidity. The product should be kept at a controlled temperature between 4 °C and 25 °C.

TRANSPORT CONDITIONS:
The product can be transported between -28 °C and 55 °C for up to 15 days without compromising its initial performance.

VALIDITY:
If stored according to the manufacturer's instructions, the product is valid for 24 months from the date of manufacture. Do not use after the expiry date.
Once opened, if stored according to the manufacturer's recommendations, the product remains stable for 24 months.
Do not use after the expiry date.

MODE OF DISPOSAL:
Do not reuse the product and dispose it in accordance with local legislation.

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PRIZMA 3D GLAZE
VERNIZ FOTOPOLIMERIZÁVEL

PRIZMA 3D GLAZE
BARNIZ FOTOPOLIMERIZABLE

PRIZMA 3D GLAZE
LIGHT-CURING VARNISH

MADE IN BRAZIL maquira Dental Group

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