PRIZMA 3D GLAZE - LIGHT-CURING VARNISH

EN

DESCRIPTION
PriZma 3D Glaze is a light-curing varnish designed for use on resins produced by 3D printing. Its application seals any porosities of the material, stabilize the surface and provide the desired shine. After application, an abrasion-resistant layer is formed, giving high shine and enhances the surface finish.

PRESENTATION
1 bottle of PriZma 3D Glaze containing 6 mL.

INDICATION FOR USE
Monomer-based varnish, indicated to sealing surface porosities of resins used in the 3D printing of complete and partial dentures. Its application promotes surface sealing, providing increased abrasion resistance and the a glossy finish to the material.

COMPOSITION
Methacrylic monomers, primers, co-initiators, stabilizers, inhibitors, fluorescent agents, and fillers.

MODE OF USE
In the dental office:
1. **Product preparation:** Properly prepare the material that will receive PriZma 3D Glaze.
2. **Product application:** Apply an even layer of PriZma 3D Glaze with a brush, using unidirectional strokes in order to avoid the incorporation of bubbles.
3. **Initial photopolymerization:** The PriZma 3D Glaze should be light-cured with blue light, using a light-curing device for 60 seconds for each layer applied.
4. **Final photopolymerization with oxygen inhibitor (glycerin):** At the end of the procedure, the application of an oxygen inhibitor (glycerin) is recommended. Immerse the characterized device in glycerin and perform the final polymerization in a UV light chamber for 3 to 5 minutes, ensuring complete polymerization.

Resin diluent for characterization:
1. **Surface preparation:** To improve the mechanical adhesion between the surface of the part and the pigment, it is advisable to blast aluminum oxide in the area to be characterized. After blasting, clean the surface with isopropyl alcohol and wait for complete evaporation. It is recommended to apply a thin layer of PriZma 3D Glaze prior to characterization with resinous pigment to promote adhesion.
2. **Pigment application:** Apply the resinous pigment with the help of a soft-bristled brush, in thin layers, until the desired aesthetic effects and coloring are achieved. If necessary, PriZma 3D Glaze can be used as a thinning agent in small quantities, just to facilitate application and allow resin shaping.
3. **Photopolymerization of the pigment:** After each layer of pigment applied, perform the photopolymerization with blue light, using a light-curing device for 40 seconds. Repeat the application and light-curing process until the complete characterization of the part is complete. After completing all the characterization and light curing of all pigment layers, perform the final polymerization in a UV light chamber for 3 to 5 minutes.
4. **Product application:** Only after the characterization is fully finalized and polymerized, apply a uniform layer of PriZma 3D Glaze with a brush, using movements in a single direction, in order to avoid the incorporation of bubbles.
5. **Initial photopolymerization of the product:** The photopolymerization of the PriZma 3D Glaze with blue light should be performed, using a light-curing device for 60 seconds for each layer applied.
6. **Final photopolymerization with oxygen inhibitor (glycerin):** At the end of the procedure, the application of an oxygen inhibitor (glycerin) is recommended. Immerse the characterized device in glycerin and perform the final polymerization in a UV light chamber for 3 to 5 minutes, ensuring complete polymerization.

WARNINGS
- Exclusive use of the dental surgeon and dental technician.
- Do not reuse empty packaging.
- Keep the product out of the reach of children.
- This material is manufactured for dental use only and must be handled in accordance with the instructions for use. The manufacturer is not responsible for damage caused by other uses or by incorrect handling not described in these instructions for use. In addition, the user is obliged to check, before use and under his responsibility, whether this material is compatible with the intended use, especially when this use is not indicated in this indications for use.
- Strictly follow the procedures provided in the instructions for use.
- Do not allow loose polymerized objects or fragments to be swallowed or aspirated.
- Contact of the uncured product with the skin, mucous membranes and eyes should be avoided. When not yet polymerized, it can cause a slightly irritating effect.
- Single-use product. After photoactivation, the product cannot be restored to its original consistency.
- Reprocessing is prohibited.
- Do not expose the material directly to sunlight.
- After use, the perfect sealing of the packaging must be observed to prevent the material from being exposed to light.
- Do not use the material if the packaging is damaged.
- PriZma 3D Glaze light-curing varnish contains methacrylates and other substances that may cause irritation or sensitization. To ensure safe handling, always use dental Personal Protective Equipment (PPE), such as gloves, mask, protective eyewear and a gown.
- Do not use the product after the expiration date.

PRECAUTIONS
- This product must be used strictly in accordance with the instructions for use.
- In case of allergic reactions, discontinue use immediately. If necessary, seek medical advice.
- In case of storage in refrigerated systems, before application, the material should be left at room temperature so that it reaches adequate fluidity.
- Each UV light chamber has different light intensity; therefore, the polymerization time must be tested beforehand. If the time initially used is not sufficient for the clinical case, it is recommended to keep the material in a UV light chamber for an additional period of time, in order to ensure complete polymerization.
- During the use of light-curing devices, always use appropriate eye protection.
- To avoid allergic reactions to methacrylate monomer, it is recommended to continuous use of protective gloves is recommended.
- If there is contact with the eyes, rinse them thoroughly with running water and seek immediate medical attention.
- Avoid contact with the oral mucosa in order to avoid the inflammatory reactions. In accidental contacts, washing with plenty of water and medical assistance is recommended.
- PriZma 3D Glaze light-curing varnish should not be ingested or aspirated. In case of accidental ingestion or aspiration, seek medical attention immediately.
- Use a light-curing device with a minimum power of 1200 mW/cm². On devices with lower power, it is recommended to increase the photoactivation time.

CONTRAINDICATIONS
PriZma 3D Glaze should not be used in patients with known hypersensitivity (allergy) to any of the components of the formulation.

ADVERSE EFFECT
To date, no adverse effects have been identified from the product PriZma 3D Glaze. However, there are reports in the literature regarding allergic reactions involving dental products based on methacrylic monomers.

STORAGE CONDITIONS
Store PriZma 3D Glaze in its original container and in a clean, dry place. Keep the material tightly closed and without exposure to light or high humidity. The product should be kept at a controlled temperature between 4 °C and 25 °C.

TRANSPORT CONDITIONS
The product must be transported in its original packaging, in a clean and dry environment. The product may be transported between -28 °C and 55 °C for up to 15 days without compromising its initial performance.

VALIDITY
If stored as instructed by the manufacturer, the product has a shelf life of 24 months from the date of manufacture.
Once opened, if stored according to the manufacturer's recommendations, the product remains stable for 24 months.
Do not use after the expiration date.

MODE OF DISPOSAL
Do not reuse the product and dispose of in accordance with local legislation.

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