PRIZMA 3D MAKEUP - LIGHT-CURING PIGMENT FOR RESIN CHARACTERIZATION **EN**

DESCRIPTION
PriZma 3D Makeup is a light-curing resin pigment, available in different colors, developed to enable the aesthetic characterization (makeup) of composite resins and 3D-printed parts. It is available in the colors Intense A, Intense B, Terracota, Honey, Vanilla, Purple, Acqua, Lunar, Red, Pink, Yellow and Black.

PRESENTATIONS
- 1 syringe of PriZma 3D Makeup containing 2 g.
- Set containing 8 PriZma 3D Makeup syringes (2 g each), 1 bottle of PriZma 3D Glaze, and 25 tips.

INDICATION FOR USE
PriZma 3D Makeup is indicated for:
- Makeup of 3D resin printed provisionals;
- Chromatic adaptation of teeth in complete dentures;
- Color correction of artificial teeth;
- Characterization of partial dentures and artificial gingiva;
- Characterization of ceramic materials or zirconia-reinforced composites, conventional composites and hybrid ceramics.

COMPOSITION
Methacrylic monomers, primers, co-initiators, preservatives, fillers, and pigments.

HOW TO USE
General Care:
Before applying PriZma 3D Makeup light-curing pigments, the surface must be properly prepared and the occlusion and contact points must be correctly adjusted.

Application:
1. **Surface preparation:** To improve the mechanical adhesion between the surface of the restoration and the pigment, it is advisable to sandblast with aluminum oxide in the area to be characterized. After sandblasting, clean the surface with isopropyl alcohol and wait for complete evaporation. It is recommended to apply PriZma 3D Glaze prior to characterization with resinous pigment to promote adhesion.
2. **Pigment application:** Apply PriZma 3D Makeup resin pigment with a soft-bristled brush, in thin layers, until you achieve the desired aesthetic effect and shade. If necessary, PriZma 3D Glaze may be used as a pigment thinning agent.
3. **Photopolymerization of the pigment:** After each layer of pigment applied, perform light-curing with blue light, using a light-curing unit for 40 seconds. Repeat the application and light-curing process until the complete characterization of the restoration is complete. After completion of all characterization and light-curing of all pigment layers, perform the final curing in a UV light chamber for 3 to 5 minutes.
4. **Application of PriZma 3D Glaze:** Only after the characterization is fully finished and photopolymerized, apply a uniform layer of product with a brush, using movements unidirectional strokes, in order to avoid the bubble formation.
5. **Photopolymerization of PriZma 3D Glaze:** Perform the photopolymerization of the product with blue light, using a light-curing unit for 60 seconds per applied layer.
6. **Final photopolymerization with oxygen inhibitor (glycerin):** At the end of the procedure, the use of an oxygen inhibitor (glycerin) is recommended. Immerse the characterized restoration in glycerin and perform the final polymerization in a UV light chamber for 3 to 5 minutes, ensuring complete polymerization.

WARNINGS
- The product is intended for professional use.
- This material has been manufactured exclusively for dental use and must be handled in accordance with the instructions for use. The manufacturer is not responsible for damage caused by other uses or misuse, unintended use not described in this manual. In addition, the user must check, before use and under their own responsibility, that this material is compatible with the intended use, especially when this use is not indicated in this indication of use.
- The professional responsible for the application must have experience in characterization of prosthetic restorations and dental medical devices, in order to ensure that the product is used as indicated and achieves the expected performance.
- Keep the material out of the reach of children.
- Do not reuse empty packaging.
- Do not use the material if the packaging is damaged.
- Do not use the product after the expiration date.
- Replace the cap on the syringe immediately after each application.
- Use a new application tip with each use.
- Single-use product. After photoactivation, the product does not return to its initial consistency.
- Reprocessing is prohibited.
- Do not expose the material directly to sunlight.
- The product should not come into contact with eugenol-containing substances, as this component inhibits the photopolymerization process and compromises the clinical performance of the pigment.
- Do not allow loose polymerized objects and parts to be swallowed or aspirated.
- Contact of the uncured product with the skin, mucous membranes and eyes should be avoided. When not yet polymerized, it can cause a slightly irritating effect.
- After use, the tight closure of the packaging must be observed to prevent the material from being exposed to light.
- PriZma 3D Makeup resin pigment contains methacrylates and other substances that can cause irritation or sensitization. To ensure safe handling, always use dental Personal Protective Equipment (PPE), such as gloves, mask, goggles and apron.
- The tips must be disposed of immediately after each use, in a container suitable for waste from health services, and their reuse is prohibited.
- Light curing must be performed properly in order to ensure the durability and color stability of the product.

PRECAUTIONS
- This product must be used strictly in accordance with the instructions for use.
- In case of allergic reactions, discontinue use immediately. If necessary, seek medical advice.
- If the product is stored in refrigerated systems, before application, it must be allowed to reach room temperature until it reaches adequate flow is achieved.
- Each UV light chamber has different light intensity; therefore, the curing time must be validated beforehand. If the time initially used is not sufficient for the clinical case, it is recommended to keep the material in a UV light chamber for an additional period of time, in order to ensure complete polymerization.
- During the use of light-curing devices, always use appropriate eye protection.
- To avoid allergic reactions to methacrylate monomer, it is recommended to continuous use of protective gloves is recommended.
- If there is contact with the eyes, rinse them thoroughly with running water and seek immediate medical attention.
- Use a light-curing unit with a minimum power of 1200 mW/cm². On devices with lower power, it is recommended to increase the photoactivation time to 60 seconds.

CONTRAINDICATION
PriZma 3D Makeup should not be used in patients with known hypersensitivity (allergy) to any component of the formulation.

ADVERSE EFFECT
To date, no adverse effects have been identified for PriZma 3D Makeup when used as intended PriZma 3D Makeup. However, the literature reports regarding allergic reactions involving dental products based on methacrylic monomers.

STORAGE CONDITIONS
Store PriZma 3D Makeup in its original container and in a clean, dry place. Keep the material tightly closed and protect the material from light or high humidity. The product should be kept at a controlled temperature between 4 °C and 25 °C.

CONDITIONS FOR CARRIAGE
The product must be transported in its original packaging, in a clean and dry environment. The product may be transported for 15 days between -28 °C and 55 °C without compromising its integrity.

SHELF LIFE
If stored according to the manufacturer's instructions by the manufacturer, the product has a shelf life of 24 months from the date of manufacture.
Once opened, if stored according to the manufacturer's recommendations, the product remains stable for 24 months.
Do not use after the expiration date.

DISPOSAL MODE
Dispose of the product in accordance with applicable local regulations.
It is recommended that the professional apply good clinical practices for the disposal of medical devices, according to current health standards.

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Light-curing pigment for resin characterization

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